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(International Society of Thrombosis and Haemostasis)

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Bibliographie

1. Parenteral Anticoagulants – Antithrombotic Therapy and Prevention of Thrombosis, 9th ed: American College of Chest Physicians Evidence-Based Clinical Practice Guidelines Chest 2012;141:e24S-e43S
2. Belk *et al.* A comparison of red blood cell transfusion utilization between anti-activated factor X and activated partial thromboplastin monitoring in patients receiving unfractionated heparin. J Thromb Haemost 2016; 14: 2148–57.
3. Vandiver *et al.* Antifactor Xa levels vs Activated Partial Thromboplastin Time for Monitoring Unfractionated Heparin. Pharmacotherapy 2012;32:546-558
4. Levine MN *et al.* A randomized trial comparing activated thromboplastin time with heparin assay in patients with acute venous thromboembolism requiring large daily doses of heparin. Arch Intern Med. 1994 Jan 10;154(1):49-56.
5. Guervil *et al.* Activated Partial Thromboplastin Time vs Antifactor Xa Heparin Assay in Monitoring Unfractionated Heparin by Continuous Intravenous Infusion. Ann Pharmacother. 2011 Jul; 45(7-8):861-8.
6. Delmas *et al.* Anticoagulation Monitoring Under ECMO Support: A Comparative Study Between the Activated Coagulation Time and the Anti-Xa Activity Assay J Int Care Med 2018; DOI: 10.1177/0885066618776937

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Diagnostica Stago S.A.S.
RCS Nanterre B305 151 409
3, allée Thérèse
92600 Asnières sur Seine
France
Tel.: +33 (0)1 46 88 20 20
Fax: +33 (0)1 47 91 08 91
webmaster@stago.com
www.stago.com